

Mediwish Surgical Instrument Tray Corner Protectors:

Manufacturer Instructions

Protecting Your Instruments and Streamlining Workflows

This document provides comprehensive instructions for using Mediwish Surgical Instrument Tray Corner Protectors (Art. STCPB40). These innovative corner protectors safeguard your valuable instruments by preventing rips, tears, and holes in sterilization wraps caused by sharp tray corners, improper handling, or heavy instrument loads.



Applications:

- Protection during sterilization processes
- Safeguarding instruments during transport and storage
- Minimizing damage from improper handling

Key Features:

- **Flexible Design:** Adapts to various instrument tray corner shapes and sizes.
- **Durable Material:** Made from high-quality, medical-grade silicone for long-lasting performance.
- **Sterilization Compatibility:** Approved for various sterilization methods, including:
 - Autoclave: 273°F (134°C) for 30 minutes
 - Ethylene Oxide
 - Ozone: 97°F (36°C) for 4.5 hours
 - Sterrad Standard Cycle: 28 minutes

Cleaning Recommendations:

- Utilize a Central Sterile Washer for optimal cleaning and disinfection.
- Alternatively, follow standard protocols for cleaning medical-grade silicone products.

For further information or technical specifications, please refer to the following resources:

- [Mediwish Brochure](#): Offers a comprehensive overview of the product's features and benefits.

Mediwish Customer Support:

- Phone: +86 (551) 65557089
- Email: info@mediwish.com