



MEDIWISH

Medical product supplier

PRODUCT CATALOGUE

•STERILIZATION PACKAGE

•MEDICAL PRODUCTS

•DENTAL PRODUCTS

Mediwish Co., LTD
Address: The Crossing of Tianshui Rd South and
Jinchi Rd, Hefei Luyang Industrial Zone
WatsApp: +8618955119211
Office phone: +86 551 65557089
E-mail: ying@mediwish.com
www.mediwish.com
www.mediwish.net



安徽美志斯工贸有限公司



Thank you career path for your help, the most sincere best wishes to you, make a beautiful wish, let all of your wishes had come true, all pray for efficacy, all promise to deliver, can you see all the joy and happiness

感事路上有的助，最真的祝福送，下美的心愿，所有的愿望都，所有的祈求都，所有的言都，所有的快和幸福都能看得

Co-Founder
Ketty Zhang

Contents

Sterilization Packaging

- Self Seal Pouches / p.3
- Sterilization reels / p.4
- Flat pouches / p.5
- Gusset pouches / p.6
- Paper gusset pouches / p.7
- Tyvek pouches & reels / p.8-9
- Ware basket / p.10
- Sterilization equipment / p.11
- Wrapping cloths / p.12-13
- Other pouches / p.14
- Chemical indicators / p.15

Seal machines

- Handle seal / p.16-17
- Rotary seal / p. 18-19

Nonwoven products

- Dental bibs / p.20-21
- Face mask / p. 22
- Disposable products / p.22-23

Dental products


- Dental instrument / p.27
- Disposable products / p.28
- Gause & dressing / p.29
- Infusion products / p.30
- Examination kit / p.31-32
- Dental sleeves / p.33
- Microbrushes / p.34
- Air water syringe tips / p.35
- Wedges, mixing tips / p.36
- Dental laboratory / p.37-38
- Dental products / p.39-40
- Instruments container / p.41
- Prophy brush & cup / p. 43
- Safety glasses / p.44
- Check-retractor / p.45
- Ortodontic / p.46
- Bandages & plasters / p.47
- Medical instruments / p.48

PRODUCTS / Sterilization Packaging

Self Seal Sterilization Pouch

Material: made of french arjo paper/medical paper +CPP/PET Film(blue/green/white)

Handle of seal: self adhesive strip sealed. Advantage: see interior items clearly through transparent film, have good sterile barrier excellent peelability, seal quickly. Indicators: color indicators for STEAM, ETO and FO. Certificate: ISO13485, CE. Usage: are used by dentists, salons, hospitals, spas, tattoo artists and many other types of medical facilities to enable the sterilization and storage of surgical instruments, medical gloves and wound care suppliers prior to use. Packaging for medical device, surgical instruments, medical consumable and so on.



Item No.	Size	Manner of packing		Outer carton size	G.W.	N.W.
SS057100	57x100mm	200pcs/box	50boxes/ctn	35x23x31cm	7.3kg	6.6kg
SS057120	57x120mm	200pcs/box	50boxes/ctn	33x29x31cm	9.3kg	8.6kg
SS070260	70x260mm	200pcs/box	10boxes/ctn	35x29x18cm	4.5kg	4.2kg
SS083160	83x160mm	200pcs/box	10boxes/ctn	35x35x27cm	10kg	9.2kg
SS090165	90x165mm	200pcs/box	10boxes/ctn	36x35x21cm	7.5kg	6.9kg
SS090260	90x260mm	200pcs/box	10boxes/ctn	35x29x21cm	5.8kg	5.4kg
SS135290	135x290mm	200pcs/box	10boxes/ctn	35x30x30cm	9.5kg	8.7kg
SS155260	155x260mm	200pcs/box	10boxes/ctn	35x29x34cm	10kg	9.3kg
SS190360	190x360mm	200pcs/box	5boxes/ctn	39x35x22cm	8.5kg	9.5kg
SS230395	230x395mm	200pcs/box	5boxes/ctn	42x35x25cm	11.2kg	10.5kg
SS300395	300x395mm	200pcs/box	5boxes/ctn	42x35x32cm	14.5kg	13.7kg
SS300430	300x430mm	200pcs/box	5boxes/ctn	46x35x32cm	16.2kg	15.1kg

PRODUCTS / Sterilization Packaging

Sterilization flat reel

Mediwish heat-sealable sterilization reels are constructed from transparent PET/PP multi-layers copolymer film and medical grade paper. Process indicators for steam and ethylene oxide sterilization are applied on the paper surface of the reel and help to differentiate between processed and unprocessed packages.

All process indicators are water based and non-toxic inks and meet or exceed the requirements of ISO 11140-1 standard and provide accurate and clear color change after sterilization process.

Wide range of standard sizes allows optimum choice of a correct-sized pack for each item and for double packing purposes and as a requirement of EN ISO 11607 all equipments, machines and processes used for Mediwish sterilization reel manufacturing, process validation has been conducted.

Product Description

- Superior barrier with 60gsm or 70gsm medical grade paper conforming EN 868-3
- Transparent, reinforced multilayer co-polymer PET/PP film, ISO 11140-1 certified water based, non-toxic and accurate process indicator
- Three independent seal lines
- Excellent gas permeability and microbial barrier properties
- Can be easily sealed with all sealing machines
- Manufactured under validated production process.



Sterilization reel with gusset

- All printing is made inside the 11 mm. sealing area to prevent pigment migration by using water based inks which are appropriate with EN 868-5,
- The product inside packaging, can be seen easily because of transparent, multi-layer co-polymer film
- High-resistant, strengthened, non-toxic film which provides aseptic presentation.
- Folding height is appropriate to standards
- High-resistant, strengthened, non-toxic film.
- High-barrier property first class medical craft paper in 60 gr/m2, 70 gr/m2
- Can be easily closed with any kind of set-top closing device
- Appropriate to the norms of non-toxic process indicators ISO 11140
- Easily opening, non-tearing with indicated opening direction
- In case of need, formaldehyde indicator can be applied.



Item No.	Size	CTN	Outer carton size	N.W.	G.W.
Code	cm	rolls	cm	kg	
GR075100	75mmx100mx25mm	16rolls	40x40x32cm	18kg	19kg
GR100100	100mmx100mx50mm	12rolls	40x40x32cm	18kg	19kg
GR150100	150mmx100mx55mm	8rolls	40x40x32cm	18kg	19kg
GR200100	200mmx100mx65mm	4rolls	40x40x22cm	12kg	13kg
GR250100	250mmx100mx70mm	4rolls	40x40x27cm	14kg	15kg
GR300100	300mmx100mx80mm	4rolls	40x40x32cm	18kg	19kg
GR350100	350mmx100mx80mm	4rolls	40x40x37cm	20kg	21kg
GR400100	400mmx100mx100mm	4rolls	40x40x42cm	23kg	24kg

PRODUCTS / Sterilization Packaging

Sterilization flat pouch

Mediwish sterilization pouches are constructed from transparent PET/PP multi-layers copolymer film and medical kraft paper. Process indicators for steam and ethylene oxide sterilization are applied on the paper surface of the reel and help to differentiate between processed and unprocessed packages.

All process indicators are water based and non-toxic inks and meet or exceed the requirements of ISO 11140-1 standard and provide accurate and clear color change after sterilization process.

Wide range of standard sizes allows optimum choice of a correct-sized pack for each item and for double packing purposes and as a requirement of EN ISO 11607 all equipment, machines and processes used for Mediwish sterilization pouch manufacturing, process validation has been conducted.

Product Description

- Superior barrier with 60gsm or 70gsm medical grade paper

- Transparent, reinforced multilayer co-polymer film

- ISO 11140-1 certified water based, non-toxic and accurate process indicator

- Three independent seal lines

- Can be easily sealed with all sealing machines



PRODUCTS / Sterilization Packaging

Sterilization Gusset pouch

It is composed of medical sterilization paper and delaminated, transparent and composite films.

The pouch being three-dimensional increases the volume, and reduces the packing cost effectively. There is a folded gusset on each side of the pouch which enables a bigger inner space than the flat pouch to load the materials to be sterilized.

So the pouch is suitable for thick objects' sterilization

The pouch provides sterilization color-changed instructions, and the sterilization condition is clear. There are two sterilization indicators on the pouch including ethylene oxide and steam, which will change colors obviously when the sterilization process is finished, and are able to indicate if the sterilization process is enough and can distinguish the articles from sterilized to non-sterilized.


Explosion-proof edges can ensure safety.

Material Structure:

cpp+pet+dialytic paper

film's color: white, blue, green, purple

Sterilization/antiseptic way: ETO, STEAM, FORM



PRODUCTS / Sterilization Packaging

Paper gusset pouch

Material: France Arjo Wiggins / Sweden Billerud (60g/m²) Medical high-temperature dialysis paper. Made by pure medical dialysis paper, environmentally friendly. Reducing the packing cost of unit product by the gusseted design. The indicator's color has significant change and no change back. The initial blue will turn into black based on steam sterilization method; and will turn brown from initial pink based on ETO sterilization method. The simple heat-sealing machine can be used to seal easily and fast. Both size and print content can be produced and designed by the customer's request.

Item No.	Size	Manner of packing	Outer carton size	N.W.	G.W.
LT090125	90X50X125	3000pcs/ctn	30x28x28cm	7kg	8kg
LT090150	90X50X150	3000pcs/ctn	33x30x28cm	8.3kg	9.3kg
LT090250	90X50X250	1500pcs/ctn	30x27x28cm	7kg	8kg
LT110190	110X30X190	1000pcs/ctn	23x21x28cm	3.5kg	4.5kg
LT040250	140X75X250	1000pcs/ctn	31x27x28cm	7kg	8kg
LT190250	190X65X250	1000pcs/ctn	41x27x28cm	8.3g	9.3kg
LT140330	140X50X330	1000pcs/ctn	35x31x28cm	6.3kg	7.3kg
LT180330	180X95X380	500pcs/ctn	40x20x28cm	6.8kg	7.8kg
LT160380	160X40X380	1000pcs/ctn	40x35x28cm	9.5kg	10.5kg
LT100150	100X50X150	3000pcs/ctn	33x33x28cm	9kg	10kg
LT190330	190X65X330	500pcs/ctn	35x21x28cm	5.5kg	6.5kg
LT110180	110X45X180	1500pcs/ctn	36x21x28cm	5.6kg	6.6kg
LT125250	125X90X250	1000pcs/ctn	28x27x28cm	7kg	8kg
LT125150	125X30X150	2000pcs/ctn	33x28x28cm	6kg	7kg



PRODUCTS / Sterilization Packaging

Tyvek sterilization roll



Tyvek is a rather unique material in the sterile packaging classification. It is known for its outstanding strength, durability and tear resistance. As well as industrial use for EO and Gamma sterilization, Tyvek is also used for H₂O₂ gas plasma sterilization at hospitals. MediWish Tyvek sterilization reels are constructed from transparent multilayer PET/PE copolymer film and uncoated Tyvek web in compliance with EN 868-9. Water based, non-toxic process indicator for hydrogen peroxide sterilization complying with ISO 11140-1 is applied on the Tyvek surface and helps to differentiate between processed and unprocessed packages.

Product Description- Made of transparent PET/PE and Tyvek; Excellent microbial barrier properties; Easy peel and aseptic presentation; Wide sealing temperature window, Easy seal with old sealing machines, Clear and accurate color change from red to blue, Suitable for plasma and ozone sterilization

Tyvek Sterilization Pouch



Constructed from 100% high density polyethylene Dupont Tyvek paper and PET/PP laminated film. Transparent multilayer laminated film allows easy identification. Dual band seal creates three independent barriers to contamination. Strong seal strength, visible sealing lines. All imprints are located outside the packing area to prevent ink migration to the products. Water based and non toxic. Compatible with all Low Temperature Hydrogen Peroxide sterilization systems. Easy to peel off.

Size range		Components		Sterilization
length	width			
10-100cm	5cm-50cm	Standard component	1073B coating TYVEK layer+ PE/PET laminated film (white)	Gamma
		Optional component	1073B un-coating TYVEK layer+ PE/PET laminated film (white)	
			1059B coating TYVEK layer+ PE/PET laminated film (white)	ETO Plasma
			1059B un-coating TYVEK layer+ PE/PET Medical laminated film (white)	

Tyvek Header Bag



Material :
Tyvek+composite film

Features :
Superior bacterial barrier ability, professional design makes this type of packaging resistance to the penetration of bacterial spores. Good permeability, Easy to use, clean, peel, avoid secondary pollution products. Suitable for ethylene oxide, gamma ray, plasma and other sterilization methods. Specification is not limited, Can be produced according to the customer's desired size.

Apply: Applicable to the higher cost, higher risk of two categories, three categories of products. Disposable medical device package, such as anesthesia bag, operation bag, etc.




PRODUCTS / Sterilization Packaging

Metalized Sterilization Pouch (Al-pouch)



Aluminum foil bag.

1. materials :PET/AL/PE material
2. size : 75x95mm or 60x80mm or others
3. thickness :100mic.
4. food grade material colorful printing, sticker and printing are available also can put zipper and tear notch...
- 1.this bag can be used in kinds of industries
- 4 . bag shape can be flat with zipper , stand up with zipper .(based on client's need)



EtO Sterilization Pouches



Heat-sealable sterilization pouches special designed for ethylene oxide (EtO) sterilization method. EtO pouches ensure reliable protection against contamination with bacteria from time of sterilization until use of the sterile medical device.

EtO sterilization pouches are constructed from transparent PET/PE multilayer copolymer film and ethylene oxide grade paper. Water based, non-toxic process indicators for ethylene oxide sterilization complying with ISO 11140-1 is applied on the paper surface and helps to differentiate between unprocessed and processed packages.

EtO sterilization pouches are available in different sizes and can be customized upon request.



Heat Sealing Paper-Paper Pouch



Sterilization Method

Steam, EO sterilization, Plasma, Gamma Ray

Apply to Widely used in hospitals, dental clinics, laboratories and medical disposable manufactures

Excellent Peelable seal and clean peel,

Excellent bacteria-barrier properties

Manufactured to individual customer specifications




PRODUCTS /Ware Baskets 304ss

Wire Baskets Ø 1 mm; 304 ss wire mesh; cell dia 4 mm



Our wire mesh and perforated sheet baskets can be manufactured in any size.
Please contact for detailed information

Endoscope baskets Ø 1 mm; 304 ss wire mesh; cell dia 4 mm





Laundry baskets



Wire baskets with Led



Accessories for cleaning / Brushes



PRODUCTS / Sterilization Equipment

Medical washer disinfectant



Washer for processing instruments in dental surgeries, Profi water softener
Wash cabinet made from high quality stainless steel, High performance circulation pump, 7 kinds of washing programs
Capacity: 65L; Wash carts can loading 16 DIN trays and dental devices; Unique double pump structure and double stainless filter make besmirch separation and automatic exclusion; Polishing agent dispenser; Only for liquid detergent
Noise: ≤50dB, Weight: ≤43kg

Ultrasonic cleaner



Capacity: 1L-10L (Inner chamber: 300×240×150mm)
Frequency: 40 kHz
Power: 240W
Temperature: up to 80°C
Adjustable Time: 1-99 minutes.
Features: PC control, auto-checking frequency, Adjustable Time, Temperature adjustable, Count down with digital display, Display current temperature and time timely.

Handpiece Lubrication Machine



Description
Handpiece Lubricator System
Three connectors
Working pressure: 0.3-0.6Mpa
Voltage: 220V/50Hz
Air flow: 60L/min
Weight: 7.2kgs
Package size:
282×300×370(mm)
165×80mm S#

Water Distiller



Specifications:
Output: Over 6.5 Gallons per Day (26 Liters)
Size: 8 in. (20 cm) width X 14 1/4 in. (36 cm) height
Weight: 7.7 lbs. (3.5 Kg)
Boiling Chamber Capacity: 1 Gal. (4 L.)
Collection Container Capacity: 1 Gal. (4 L.)
Power Source: 220-240V/60Hz
Power Consumption: 750W (total)

Sterilization tray



sterilization tray
Specification:
Shallow type:
204×130mm XL#
190×116mm L#
178×104mm S#
Deep type:
243×125mm XL#
203×102mm L#
165×80mm S#

Cleaning perforate tray



Dental Dry heat autoclave



Class B autoclave with European standard EN13060.
.High precision control system for more perfect sterilization.
.Friendly and intuitive user interface with LCD screen, displays the information about cycle process by a large graphic all the time. Have 4 display languages for the user choosing, which is English, German, Spanish and Chinese. A built-in independent steam generator, making the sterilization fast and completely.
Model-MZS-18L-C
Voltage(V/Hz) 110V-240V/50Hz-60Hz; Power(W) 1800W
Chamber size (mm) 247×350; Packing size (mm) 690×590×530
Net(kg)-48; Gross(kg)-53



PRODUCTS / Sterilization Packaging

Crepe Wrapping Paper



Item No.	Size	Manner of packing	Outer carton size	N.W.	G.W.
CP030030	30X30	1000pcs/ctn	32X32X18cm	5.4kg	6.kg
CP040040	40X40	1000pcs/ctn	42X42X18cm	9.6kg	10.1kg
CP045045	45X45	1000pcs/ctn	47X47X18cm	12.2kg	12.7kg
CP050050	50X50	800pcs/ctn	52X52X14cm	15kg	15.5kg
CP060060	60X60	500pcs/ctn	62X36X17cm	11kg	11.5kg
CP090090	90X90	500pcs/ctn	92X53X17cm	24.3kg	24.8kg
CP100100	100X100	300pcs/ctn	102X57X11cm	18kg	18.5kg
CP120120	120X120	300pcs/ctn	122X67X11cm	26kg	26.5kg

Wide range of wraps for different purposes.

crepe wrapping paper are made of 100% cellulose fibres and can be used as either inner or outer wrapping. Crepe wraps are suitable for steam, ethylene oxide, formaldehyde and radiation sterilizations methods and are a reliable solution for preventing cross contamination with bacteria.

Wrapping papers offered by Mediwish are in compliance with all international standards.

PRODUCTS / Sterilization Packaging

MEDICAL NON-WOVEN WRAPPING CLOTH

SMS, SMMS and SMMMS non-woven fabric are an extraordinarily water proof and resilient non-woven material made of polypropylene fabric layers. They can be used either as inner or outer wrapping for medical devices and trays and available in 40g regular weight, 50g heavy duty and 60g multi-pac.

SMS, SMMS and SMMMS wraps as single or double wrapped are suitable for steam, ethylene oxide, hydrogen peroxide plasma and irradiation sterilization methods. The memory-free feature before and after sterilization, ensures good aseptic presentation.



PRODUCTS / Sterilization Packaging

Easy tear pouch



Protective tubes

Some medical and dental devices which have sharp end when in sterilization process usually break the pouches or wrapping cloth the protective sleeve protect the sterilization package in case of breaking package



Tray for sterilization

Medical paper tray

1. medical PE coated paper
 2. usage for hospital, dental clinic
 3. EO, STEAM and GAMMA
- Sizes: 205x110x10mm, 225x140x-15mm, 255x135x15mm



Mid-sealing poly pouch

Material: plastic or other lamintion
Size & logo: based on client's requirement
High temperature resistance
 The pouch could be made of different material structures, suitable for packing different kinds of products.



Kraft-Paper Self Sealing Pouch

Paper paper flat pouch, made from medical paper with glue seal, seld seal type for steam, eo, hot air sterilization process



Steam indicator label

ETO STICKER
 18X25MM
 50PCS/SHEET
 155x200MM



PRODUCTS / Chemical indicators

Class 1 Steam Indicator



EO gas indicator



Apply to monitoring EO Sterilization.
Storage: protect from light, store at room temperature (15-30degree) and away from corrosive gas.
Notice: Only used in chemical monitoring of EO sterilization, can not be in biological monitoring. Standard: Complies with ISO11140-1 EN867/CLASS 4.

Form gas indicator 1, 4 class



Bowie-Dick Test Pack



Bowie Dick Test Pack is a single use test designed to detect the presence of residual air in pre-vacuum steam sterilizers operating at 134°C for 3.5 minutes. The indicator sheet Bowie Dick Test Pack will demonstrate a uniform color change from blue to dark brown/black when proper sterilization conditions are met and no air is present.

Autoclave tape ETO tape Form tape




1.2 mmx50m
1.9 mmx50m
2.2 mmx50m
2.5 mmx50m
or custom-made size

STEAM 1, 4, 5, 6 class



can be with adhesive tapes at the back

Irradiation indicator (X-ray)




Labeling machine



Label Guns provide an easy handling and printing solution for indicator labels used for recording purpose. Indicator labels are designed to create a direct link between the sterilization process and the patient. The label gun can be used to print important traceability data on the label such as sterilization/expiry date, operator code and autoclave/cycle no on the label.

Class 1 H₂O₂ Indicator





PRODUCTS / Sealing&Cutting Devices

Sealing Machines

We manufacture medical sealing and printing machine, paper-plastic-film cutter, worktable, label printer as well as some matching auxiliaries which are used in sterilizing centre, operation sterilizing room and Dental sterilizing department, including management software, paper and plastic bag, seal test paper, seal strength tester, various type of label paper with sterilizing instructions and small scale washing equipment, etc. We manufacture over twenty different types of seal machines to meet various consumer requirements. There are two type of operational systems basically, one is manual seal another is automatic seal, both of them have different control modes from temperature controller, liquid crystal screen, touchable screen, computer, scanner, etc.

PRODUCTS / Sealing&Cutting Devices

Hand sealer MZS300

Hand sealer MZS400

I. Summary

MZS300 (400) that we produced is a manual one which is for seal bags. This product using some advanced technology such as dual heaters temperature adjustable ultra high temperature automatic protection low temperature alarm over-current protection with the advantage of fast heating rate stable performance aesthetic appearance using security easy to operate low repair rate. This mode is very popular in the hospital and clinic.

II. Technology parameters

Pressure AC220±10 50Hz

Power 500W

III. Handle Instructions

Open the box and check the machine. Be sure that the machine is intact the screw tightening. Repeated slide the cutter and wrench making sure there is no internal damage.

Put the paper-plastic bags on the bracket and then import the paper shroud.

Turn on the power, adjust the knob to a fittable temperature (we recommend 170 users can adjust the temperature according to the local climate). The indicate light lit, when the light is off the machine can start work.

Pull out the bags pull the wrench on the right side for 2-3s and then slide the cutter. Push the wrench and then take the bags out.

IV. Shape and weight

Size 450×390×240(mm)

Weight 9kg

V. Storage environment

Temperature -10 ~55

Humidity ≤90% R.H

Atmospheric Pressure 50kpa~110kpa

Painted or stainless steel, seal size 300 mm or 400 mm



PRODUCTS / Sealing&Cutting Devices


Rotary sealer MZSR-1

4.3" (95×55mm) colored liquid crystal screen, light-touch key-board, English operational interface, internally installed clock;

Temperature controlled by a microcomputer, accuracy $\pm 1\%$, working temperature range 60~220 ; High rate increase of temperature, 40 seconds required from room temperature to 180 ; Automatic temperature reduction performance by an air removal system which is connected to a temperature control system to reduce the waiting time from high sealing temperature to lower sealing temperature. Adjustable fixed-force system, suitable for sealing paper-plastic bags, 3D paper-plastic bags and paper-paper bags; Advanced flat ceramic heating components, high-temperature stability, long life expectancy and high heat efficiency.

Automatic fault alarm indicator, intelligent auto-detection will alarm and remind punctually when a fault appears. Technical parameters Sealing speed 10m/min; automatic seal printing test using light-control technique. Seal width 12mm; Sealing strength meets YY/T 0698.5-2009 requirement Seal margin 0~35mm adjustable; machine shell material: spray carbon steel or 304 stainless steel Voltage: 220V 50Hz max. power: 500W

Equipment dimension: 560×260×220 mm weight: 18kg Features of Printing System: MZSR-1 can be used to print Chinese, English, numbers, and special symbols which are in "YY 0466-2003 symbols used for medical equipment labels, signs and information". It meets the requirements from health department such as, the validity of expire date, lot number, operator, etc.





Rotary sealer MZSR-2

Performance Features 4.3" (95×55mm) colored liquid crystal screen, light-touch keyboard, humanized Chinese operational interface, internally installed clock; Temperature controlled by microcomputer, accuracy $\pm 1\%$, working temperature range 60~220; High rate increase of temperature, it only needs 40 seconds to get temperature from room temperature to 180; Automatic temperature reduction performance operated by an air removing system which is connected to temperature control system to reduce the waiting time from high sealing temperature to lower sealing temperature. Adjustable fixed-force system, suitable for sealing paper-plastic bags, 3D paper-plastic bags and paper-paper bags; Automatic energy-saving during waiting time, waiting time is adjustable, intelligent recovery from waiting; restore working temperature within short time, with prolonged life time and energy saving performance. Advanced flat ceramic heating components, high-temperature stability, long life expectancy and high heat efficiency. Suitable to be connected with surrounding equipments such as multi-functional stainless paper-cutting worktable, axletree transported worktable, etc.

Technical parameters

Sealing speed 10m/min; Voltage: 220V 50Hz. Seal width 12mm; power: 500W Seal strength meets YY/T 0698.5-2009 requirement; Dimension: 490×260×136 mm. Seal margin: 0~35mm adjustable; Weight: 16KG



PRODUCTS / Sealing&Cutting Devices

Hand sealer MZS320



Sealer by hand , and with a cutting together.
Machine sealing length 300mm. Seal width 1cm.
The heating ceramic pole is from ITALY ZOPPAS.
Constant temperature 220 ± 10 °C

NET/GROSS	COLOR
7.5/7kg	painting

Hand sealer MZS-WALL



Wall type,
can hang on wall, the heating ceramic line is from ITALY ZOPPAS,
Constant temperature 220 ± 10 °C

NET/GROSS	COLOR
6.5/6kg	painting

Hand sealer MSZ100



sealer by hand , and with a cutting together.
Machine sealing length 250mm . width 1cm.

NET/GROSS	COLOR
6.5/7kg	stainless steel


Hand sealer MZS80



sealer by hand , and with a cutting together. Machine sealing length 250mm . width 1cm.

NET/GROSS	COLOR
6.5/6kg	painting white mixed with Stainless steel


Double cutter MZSC7802



Double layer paper cutter
Size (L×W×H) 780×340×380

NET/GROSS	COLOR
4.5/5kg	Stainless steel


Single cutter MZSC7801



cutting reel pouch before rotary sealing side
Size (L×W×H) 780×340×140 mm cutting width: 750mm

NET/GROSS	COLOR
7.5/7kg	stainless steel

Rolling tray MZST



Rolling tray :rolling for taking goods.
often under with rotary machine
Size (L×W×H) 745×286×52 mm

NET/GROSS	COLOR
5.5/5kg	Stainless steel

Daily checking seal test paper



Sealing test paper was "confirmation request forming, sealing and assembly processes," the need for sealing devices running Identification (OQ) requirements according to ISO 11607-2:2006 Part 2
-Sealing pressure test whether passage or opening lines appear;
-Sealing test whether the pressure is too high or too low;
-Sealing test whether the temperature is too high or too low; -Test whether continuous or sealing leaks occur;

PRODUCTS / Dental bib

Dental bibs

Principle: using positive suction, opposite the water-proof, better to protect the user's clothing.

Specifications:

1. Waterproof, absorbent feature

2. good quality

Description: Material: 1 layer membrane + 2 layer tissue ,
1 layer membrane + 1 layer tissue ,

Size (width X length) : normal size: 13 "X 18", namely 33 X 45 cm;

Unconventional size: 33 cm x 46 cm; 33 cm x 48 cm

Packing materials: the inner packing: pressure packing into the bag. out side packing: corrugated paper.

Color: white, blue, green, purple, yellow and so on.

Packaging Details: 500 pcs/box

(125pcs / PE bags, four packages in each box)

Delivery Detail: 15-25days

Dental clip



Dispenser



PRODUCTS / Dental bib with tie

Dental bib with tie

Details of specification:

Material- 2-ply paper + 1-ply poly or 1-ply paper+ 1-ply poly

Color- White, yellow, pink, green, blue, orange, purple, silver, grey, rose, beige, pich, aqua.

Size - 44x36.5cm, Weight Paper: 16-19gsm, Film: 12-14gsm. Packaging 600pcs/bag, 2400pcs/ctn. Storage

Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse. Sterilization - Sterilized by ethylene



With printing

Material	Color	size	Packaging	certification	remark
2-ply paper + 1-ply poly	blue, white, green, yellow, lavender, pink	44x36.5cm	600pcs/bag, 2400pcs/ctn	CE0197, ISO13485, ISO9001	specifications can be according to customers' requirements.

PRODUCTS / Face mask

Non-woven face mask, 3 ply, ear-loop



Specifications:
PFE> 99%@0.1 micron
Resistance to penetration by synthetic blood—160mmHg
3-ply, 3-folder comfortable design. Fluid resistant perforated film (PPF) as outer layer. Meltblown filtration middle layer. Silk-like, ultra-soft cellulosic white inner layer. Latex free, fiberglass free hypoallergenic material. Full length aluminum nose piece.

Non-woven face mask, 3 ply, with tie



Specifications:
PP non woven surgical face mask 3 plies with ties
Color: blue, green, white, pink etc
Any size available
50pcs/inbox, 2000pcs/ctn
ASTM F2100 Level 3 Barrier Procedure Face Mask
50pcs/box, 8 boxes/case
36.3*21*19.5cm

Earloop procedure mask with plastic eye shield



These pleated earloop procedure masks feature a plastic wraparound shield that helps protect the face and eyes from body fluids and pathogens for 99% average bacterial efficiency. Latex-free
With Ear loop/ Tie-on
Any size available
Packaging Details: 5pcs/bag, 10bags/cases

Non-woven activated-carbon face mask, 4 ply, ear-loop



Material-1. PP/S (100% Polypropylene Spunbond Nonwoven Fabric)+2. SMS (Polypropylene Spunbond Nonwoven Fabric + Melt-blown Nonwoven fabric + Polypropylene Spunbond Nonwoven Fabric) +3. Active carbon filter in the middle layer.
Any size available/With Ear loop/ Tie-on

Non-woven face mask with print, 3 ply, ear-loop



Disposable face mask printing details
Disposable face mask with full printing (printing all over the face mask)
1.Regular Packing: 50pcs/box, 2000pcs/ctn
2.Individual Packing: 1pc/pe bag (logo printing is available) 50pcs/box, 2000pcs/Ctn
3. With Ear loop/ Tie-on
Any size available

Disposable non-woven cone mask Respirator



Soft, flexible nose piece. Splash repellent barrier. Flared edge prevents irritation. Adjustable headband and inside soft foam cushion to adapt to various face-shapes with comfortable sealing. The high-efficient and strong separate harmful materials such as glass fiber, asbestos, etc.
Packaging Details: 5pcs/bag, 10bags/cases

PRODUCTS / Apparel & gloves

Nonwoven shoe cover



Material-PP, Great elasticity, More durable, Anti-slippery, texturized
Elastic band and contoured styling
Breathable, hygienic, Comfortable
Widely used in clinic, hospital, restaurant, food processing, beauty salon, electronic industry, etc.
Size-16*40cm, 17*38cm, 17*41cm, 17*45cm, etc, Density-30-45gr/cm²

Medical nonwoven cap

Material : PP, PE, nonwoven fabric.
Soft, breathable, comfortable, non-toxic.
It can be extend and shrink, easy fold
style offers you cost-saving, economical, elegant and good looking.
Color: rainbow color
Size: 24 inches : Weight: 12gsm
Packaging Details: 100pcs/polybag
10bags/ctn or 20bags/ctn



Isolation gown

Material: PP, SMS, MP, Tyvek
Colour White, Blue, Green
Size: S-5XL, Weight: 22-55g/m²
Packing: 5pcs/polybag, 10polybags/ctn.
The style, color, size, weight and packing can be customized



Surgical gown



Available in Standard and High Performance Surgical Gowns
Sterilized by Ethylene Oxide, for single use only
Made of non woven material for fluid resistance
Soft and comfortable disposable surgical gowns
Folded ready to use, Each gown is packed in one sealed pouch

Lab coat



Material: PP, SMS, MP, Tyvek
Colour : White, Blue, Green
Size: S-5XL
Weight: 22-55g/m²
Button/velcro/zipper front
Packing: 5pcs/polybag, 10polybags/ctn.
The style, color, size, weight and packing can be customized

Examination glove



Nitrile Gloves



PE gloves





PRODUCTS / Plastic Shoe cover



PP CPE Disposable Shoe Cover

This material of non-woven shoe cover has high degree of strength and elongation, good property of air though, water repellent, harmless and antibacterial. Non-woven Shoe Cover has been widely used in industry, agriculture, hospital and daily use.

1. Non woven material , with high intensity, non-skid, anti-tear and extensible.
2. Universal fittings, heat sealing, safe, hygienic and keeps your work place dust free.
3. Recycling , eco-friendly, comfortable, no stimulus to skin.
4. Disposable , waterproof.
5. Customised thickness, weight, color, size and styles is available.
6. High quality, lower price and best services.

Description	Specification	Weight	Pcs/Ctn	Carton Size-cm
CPE/PE Shoe Cover	15*39cm(41)	1.5g	6000	56*18*47
CPE/PE Shoe Cover	15*39cm(41)	2.0g	5000	54*18*52
CPE/PE Shoe Cover	15*39cm(41)	2.5g	4000	52*19*52
CPE/PE Shoe Cover	15*39cm(41)	3.0g	2000	55*20*38
PE Shoe Cover	15*39cm(41)	3.5g	2000	42*40*37

PRODUCTS / Non-woven products

Bed sheet disposable sterile/non-sterile



Material -1. PP/S (100% Polypropylene Spunbond Nonwoven Fabric)

2.SMS (Polypropylene Spunbond Nonwoven Fabric + Meltblown Nonwoven fabric + Polypropylene Spunbond Nonwoven Fabric)

3.PP/S+PE Film. Weight- 14gsm, 30gsm, 40gsm . 8~240gsm material

Characteristic: Durable, eco-friendly, non-toxic, breathable and flexible, Anti-alcohol, Anti-Blood, Anti-Oil, Waterproof, Acid proof, alkali proof. Size- 80*180cm , 100*180cm , 150*220cm , 200*240cm , OEM/ODM Welcomed

Color- Sky Blue, white, green, purple, or any other customized colors

Packaging- 10 pcs in PE bag, 100 pcs/box.

Storage- Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse, This product is sterilized by ethylene oxide

OEM- Material or other specifications can be according to customers' requirements. Customized Logo/brand printed. Customized packaging available.

Bed sheet in roll with perforation



Application: Disposable hospital bed roll is used for Hospital & clinic exam table and massage bed, it has perforation line for each tearing for single use

Material: Laminated paper 1-2 lay tissue with PE film ,or Nonwoven (S, SS, SMS) as per customer request. Size: minimum width 50cm x length 40m , perforation at 50cm, maximum width 110cm x length 40m , perforation at 2 meter, Or customized width ,length and perforation size. Color: blue, green, white , pink, etc.

Features: The unique technique used to bond the layers of tissue and ply securely together, which eliminates the separation from each layers. The side of the paper is absorbent, the other side of polythene is waterproof. Exam couch paper roll is a good protection for leakage when examination and massage. Remarks: Other Materials of PP Nonwoven, PE film, ONLY Tissue and nonwoven laminated polythene also can be manufactured, please contact us with material, size, inner core, quantity and package for accurate price .

Medical roll holder



Universal Paper Roll Holder

This simple east to install design

Adjustable to accommodate different size paper rolls up to 21 inches (53 cm) wide

Plugs into the table's outlet holes.

Disposable Pillow cover



Material: PP nonwoven/SPP/SMS/ PE+PP,etc. Color: white, blue, green, yellow, pink, etc. Weight 20-55g/m², Size: 50*60cm, 60*60cm, 45*75cm. Pattern plain. Features -soft, light, comfortable, breathable,dustproof, eco-friendly,etc. Usage-hospital, lab, industry, hotel, home, beauty, SPA, medical etc

PRODUCTS / Non-woven products

Dental chair cover



These full chair sleeves are made of clear plastic or SMS nonwoven and designed to easily fit for dental chairs. They provide protection against cross-contamination and minimize time spent on cleaning. Durable, convenient to use and easy to dispose of when finished. Application: Hospital, dental clinic, Medical, Cleanroom

Non-woven stretcher cover



Stretcher Cover is made from non-woven material or CPE material with elastic ring. Used to cover the Stretcher, examination couch in hospital, high-strengthened and some types are strong enough to resist a human's weight as a transfer sheet. Available in SMS NON-WOVEN, CPE or PP+PE. Available in white or blue color. Available in different size and package

Disposable head rest cover



10x10" or 10x13" Non-woven +1 layer PE 1 layer tissue+1layer PE or 2 layer tissue+1layer PE. Can be made with different texture



Poly apron



Product- PE Apron Material -Polyethylene Size : 110,70x117,80x125cm or as customer request Weight: 5-50g/pc Thickness - 0.008-0.045mm Surface - Smooth, embossed Color-Clear, white, blue, red Package-100pcs/bag, 10bags/ctn, 100pcs/box, 10boxes/ctn, 1pc/bag, 1000bags/box

PVA water soluble laundry bag



Medical PVA water soluble laundry bag Dissolvable Seam Laundry Bags (often referred to as Alginate Bags, Soluble Seam Bags or Dissolvo Bags) biodegradable, hot water soluble, dust free, Eco-friendly PROUDLY water soluble laundry bags are a cost effective solution which substantially reduces the risk of contamination and cross-infection to staff, patients and laundry personnel. Contaminated linen is placed in the bag which can then be sealed with the attached cold water soluble tie tape. The filled bag is then transferred to the laundry and out unopened into the washer. The laundry bags dissolve completely in the wash, leaving no potentially infected plastic waste. The contents of the bags need not be handled until the wash and dry cycle is complete. feature:

- 1.Avoid direct contact between the contaminated content inside the bag and the operator until the washing and drying cycle is completed, minimizing cross-infection and germs spread;
- 2.Avoid contaminated medical waste to be generated as the disposal bags completely dissolve in water, eliminating the traditional treatment by incineration or burying and reducing the treatment cost;
- 3.Greatly improve the hygiene of working environment, reduce working intension and enhance the labor protection.

Closestool cover



Size: 38*45cm
1 layer tissue+1 layer PE
or
2 layer tissue+1 layer PE

Nonwoven body washing glove Surgical scrub brush & nail cleaner



Suitable for clinic hospital patient washing body
1) Size : 15x23cm, 15x22cm, 16x23cm, 16x24cm, 17.5x25cm
2) Material : needle punched nonwoven/spunlace nonwoven/ airlaid paper laminated PE/line opening. Single sawing, over-lock sawing, glued, ultrasonic. Inner packing :100pcs/pack, 10pack/ctn Standard export carton



Medical hand scrub brush with sponge with Chlorhexidine Size: 8*5*2.1CM Sterile: Gamma Color:yellow with 4% chlorhexidine gluconate solution.

Packing: 1pc/blister, 300pcs/, 42*40*40cm, 11kg/12kg

PRODUCTS / Dental instruments

Evacuator Suction Tip (High Volume Evacuator saliva ejector), vented type



Vented
panton color
Size: 148*11mm
25pcs/bag, 50 bag/ctn

Surgical aspirator Tip



Surgical aspirating tip oral surge for dental surgery
Colour: Green
Size: 1/4", 1/8", 1/16", 3/8", types; vented/nonvented
Material: ABS
Usage: dental disposable

Surgical tip-Lamellirostral



Material-PP
packing specifications:50pc/bag, 20 bags/carton
Carton size: 38*31*27cm
Weight 4.6kgs
Color Pantone color

Evacuator Suction Tip (High Volume Evacuator saliva ejector), non-vented type



Non-Vented
panton color
Size: 148*11mm
25pcs/bag, 50 bag/ctn

Flex bend saliva ejector



bendable. 100 pcs/bag, 250 bags/ctn

Disposable Saliva Ejector



Material-Transparent PVC tube,copper interior wiring with chrome plated
Size-150*6.5mm
Packing 100pcs/bg,20bgs/ctn
Material or other specifications can be according to customers' requirements.
Customized Logo/brand printed.
Various colours available.

Disposable surgical aspirator



Articulating Paper



Specification
1. Straight or Horseshoe Shape
2. 100u/50u/30u
3. Red and blue color
4. Straight :10sheets/book, 20books/box ; Horseshoe: 6 books*20sheets
5. Thin or Thick



Dental plastic tray



Product Name	Customized Plastic Tray for Medical
Material	PMMA(Acrylic), ABS, PS, PP, PE PVC, PET, PETG, and HIPS ETC
Size	by customer demand
Color	Any color is available
Shape	Any color or Special Shape
Logo	Sticker/Screen Print etc
Packing	customized packing
Design Service	OEM/ODM Service

PRODUCTS / Dental products

Cotton rolls dispenser



112*48*120, Autoclavable, can be heated to 134 streamline design, all clear material, elegant appearance. Foam case for the entrance of the cotton rolls. Packing: 1pcs/box, 30boxes/ctn; 440*290*400mm Color: Red, Green, Blue, Black.

Denture box



Usage: Hold teeth, impression tray and so on.
Size (mm): 88x75x60, 95x75x59, 95x75x38, etc
Material: ABS, PC, PE, PP

Dental cotton roll



Specification:
8x38MM,(0.38g/pc)
10x38mm(0.45g/pc)
12x38mm(0.52g/pc)
100% absorbent cotton

It's made of the high quality 100% cotton.

Packing: 50pcs/bundle, 40bundles/bag, 8000pcs/carton.

Dental cup plastic



Sponge for cleaning



EVA EPE PE PU EPDM PVA 14cm stick 15cm total length sponge head 2x2.5cm package 5pcs/bag
2*10x50 Size by custom-made

PRODUCTS / Gause & dressing

100% Cotton gauze sponges



Cotton gauze made from pure virgin white cotton, with 100% cotton fiber filling. They provide the maximum rate and volume of absorbency combined with increased patient comfort

100% Gauze, Highly absorbent, Mesh:17*7,19*9,19*15,25*17,

Non-sterile/Sterile, CE0197,ISO13485,ISO9001 , Absorbency =3s, whiteness =82% A

Low Lint, soft and 100% all natural

Cotton filler is machine distributed evenly to minimize the bioburden


Material	Size			Packaging	OEM		
100% Cotton, Low Lint, soft and 100% all natural	2"X2" (5X5cm) 8-ply	3"X3" (7.5X7.5cm) 8-ply	4"X4" (10X10 cm) 8-ply	Package: 200 pcs/sleeve, 25 sleeve/case	Can be according to customers' requirements.	Customized Logo/brand printed.	Customized packaging available.

Disposable Sterile Gauze Ball



1. 100% Cotton, elastic ring
2. Keep all edges inside to avoid loose thread
3. Material: 100% cotton.
4. Color: white.
5. Diameter: 10mm, 15mm, 20mm, 30mm, 40mm...
6. With or without X-ray detectable thread.
7. Sterilized or non-sterile.
8. with or without elastic ring.

Cotton swab stick



- single sided medical wood stick cotton swab/cotton bud
1. White double-end cotton bud
 2. 100% natural cotton
 3. sterile/non sterile
- wood stick 7.5cm(L) Dia2.2mm
 - wood stick 15cm(L) Dia2.2mm
 - Plastic stick 7.5cm(L) Dia2.2mm
 - Plastic stick 15cm(L) Dia2.2mm
- 100pcs/bag 10bag/box
10boxs/ctn,

Alcohol swab



Specifications

1. 100% pure cotton
2. Cleanroom products.

Alcohol prep pad



1. Contents: 70% isopropyl alcohol (75% alcohol)
 2. Packet size: 50x50cm
- The common size of Alcohol Pad is 3x6 cm. But we are capable to provide various sizes of pads, such as 5 x 5 cm, 5 x 6 cm, 5 x 9 cm, and 6 x 8 cm aluminum foil pack, 3 x 6 cm, 6 x 6 cm, 12 x 20 cm non-woven or wet strength paper material.

PRODUCTS / Infusion products

1.7-2.0cc anesthesia dental cartridge syringe



Size-1.7cc-2.0cc, Color-blue
Barrel: Made of PP, Putter: Made of high grade medical ABS resin. Luer lock: Needle screwed onto the tip of syringes, rest assures more secure and the needle can be easily slipped out of the sheath, to avoid the difficulty of hand-pulling. Needle: Metric and Imperial dental needle can be applied. Sterilization: Sterilized by E.O. gas, non toxic, non pyrogenic, rest assured to use. Package: Individual blister packing, OEM is available. Application: Used in dentistry or stomatology, articaine anesthesia injection.

Disposable Dental Needle



1. Made of high quality stainless steel
 2. With large ID and high flow
 3. Sterilization: EO.
 4. Size distinguished by the needle-hub's color for clear recognition
 5. The Dental needle is available in short and long sizes for conduction and infiltration anaesthesia. It is composed of upper and lower protective cap, needle and needle hub.
- The cannula is precisely processed and siliconized to reduce pain and skin trauma for patients. The needle hubs are color coded by gauge for easier identification.
6. Procurement of all kind of special needles required by our customers

Dental Irrigation Needle



1. Needle: Made of high quality medical stainless steel.
2. Needle Tip: Blunt end, soft needle .straight or pre-bent tips resists crimping.
3. Hub: Made of medical high grade PP, universal luer lock style hub.
4. Color: According to hub's color identify needle's size, easy to distinguish and use.
5. Application: Used in dentistry or stomatology flushing and cleaning.
6. Sterilization: Sterilized by E.O. gas, non toxic, non pyrogenic, rest assured to use. carton.

Disposable Syringes



Quick Details:
Type: Luer lock, luer slip
Size: 1 ml, 3ml, 5ml, 10ml, 20ml, 50ml
Color: clear with clear plugger, clear with color plugger
Certificates: CE0197, ISO9001, ISO13485,

Safety Syringes



Properties: Injection & Puncture Instrument
Product: Safety Syringe
Component: Three parts, rubber piston, needle, barrel.
Material: Made of medical grade pp
Certificate: CE, ISO13485
Package: PE/Blister package
Brand: OEM is available
Function: After injection, the needle will retract into barrel
Color: White or Transparent
Needle: Luer lock
Size: 1ml, 2ml, 3ml, 5ml, 10ml

Sterile Infusion Set



Specifications
Disposable Sterile Transfusion Set
1. CE & ISO certified
2. High quality
3. Individual Blister pack or PE bag
4. OEM

Biohazard garbage bin bag Sharp box waste container



Material HDPE, LDPE
Size-24"x33", 26"x36", 40"x48"
50"x35" or customized, regular duty: 8-10Gallon, 12-16 Gallon, 20-30gallon, 44 gallon. Heavy duty: 8-10 Gallon, 20-30 Gallon, 44Gallon. Color-RED, BLUE, YELLOW, etc. Printing-with international biohazard symbol, or custom print. Feature-High tear and puncture resistance. seal-Flat seal or star seal packing, roll pack or loose pack. For hospital, medical clinics, and any medical



of sharps and medical waste. They can collect, store, transport and dispose of medical waste in hospital, clinic and family care.

Material- Polypropylene(PP), Paper. Volume Range 0.2L to 23L
Available Colour-Red or Yellow or Custom
The Various sharps container are designed to permit safe disposal

PRODUCTS / Dental examination kit

Dental examination kit



1. Application -used for the diagnosis and treatment of teeth, check the bad teeth in hospital and clinic.
2. Contents included: dental mirror, dental probe, forceps(tweezers), dental bib, plastic tray, saliva ejector, a pair of disposable sterile gloves, four pieces of cotton roll, a piece of tissue, cap, surgical face mask, cap for nurse.
3. Sterilized-It's sterilized by E.O. Gas and safety to be used. For single use only. Do not use if packages opened or damaged.
4. Package-The package is PE bag,the sterile gloves is packed by PE bag separately.
5. Storage-Product should be placed no more than 80% relative humidity, clean, dry well-ventilated, non-corrosive gases indoors.
6. Packaging Details: Plastic bag ,220 set per carton; Size of carton: 62*41*41cm
260 Cartons for 20' ; 620 Cartons for 40' ; G.W./N.W.:17/15KGS

Avialable compoments can include:



Components	Materials	Color
Single (or dowble)-ended Probe	ABS+Stainless steel / Carbon steel	White, green, yellow, blue, etc
Dental mirror	ABS+Glass	White, yellow, red, green, etc
Double-ended forceps (Tweezers)	ABS+Carbonsteel	White
Plastic tray	Plastic	White
Cottonroll	Cotton	White
Tissue	Paper	White
Sterile gloves	Latex	white, blue
Dental bib	One side paper, other side plastic	Blue, white
Saliva ejector	PVC	Transparent, white, yellow, ect
Cap for nurse	Non-woven	Blue, yellow, green, white, etc
Face mask	Non-woven	Blue, yellow, white, etc

Dental mirror

one-sided



double-sided



1. Used to inspect the position of bad-tooth in hospital and dental clinic.
2. The end of dental mirror is a dental spatula used to grind and mix up materials such as cement and glass ion.
- 3.Color: White, blue, yellow, green, pink. It is optional as clients' request.
4. Sterilized by EO gas and safety to be used.
5. The surface of the mirror is smooth and clean, anti-fog,no scraping,no mechanical nick, no crack and chip in the dege.The mirror is connected firmly with the mirror support,and it is level and straight.
6. Size: 160(±5)× 5(±0.5) mm; material: ABS and glass.
7. Packing material: plastic bag or blister bag

Dental probe

one-sided / double-sided



1. Used to inspect the position of bad-tooth in hospital and dental clinic.
- 3.Color: White, blue, yellow, green, pink. It is optional as clients' request.
4. Sterilized by EO gas and safety to be used.
5. Packing 100pcs/bag, 25bags/ctn
7. Packing material: plastic bag or blister bag

Dental twizers

ABC, ABC&steel, totally steel



1. Used to inspect the position of bad-tooth in hospital and dental clinic.
3. Color: White, blue, yellow, green, pink. It is optional as clients' request.
4. Sterilized by EO gas and safety to be used.
5. Material: ABS, ABC&steel, steel
7. Packing material: plastic bag or blister bag

PRODUCTS / Dental examination kits

Dental kit 2 in 1



Content:
1. Dental mirror-1pc
2. Dental probe-1pc

Material or other specifications can be according to customers' requirements.
Customized Logo/brand printed.
Various colours available.

Dental kit 3 in 1



Content:
1. Dental mirror -1pc
2. Dental probe -1pc
3. Dental tweezer -1pc

Material or other specifications can be according to customers' requirements.
Customized Logo/brand printed.
Various colours available.

Dental kit 4 in 1



Content:
1. Dental mirror -1pc
2. Dental probe -1pc
3. Dental tweezer -1pc
4. Tray -1pc

Material or other specifications can be according to customers' requirements.
Customized Logo/brand printed.

Dental kit 5 in 1



Content:
1. Dental mirror -1pc
2. Dental probe -1pc
3. Dental tweezer -1pc
4. Dental bib -1pc
5. Tray -1pc

Material or other specifications can be according to customers' requirements.
Customized Logo/brand printed.
Various colours available.

Dental set 6 in 1



Content:
1. Dental mirror -1pc
2. Dental probe -1pc
3. Dental tweezer-1pc
4. Dental bib -1pc
5. Dental gloves -1pair
6. Tray -1pc
Material or other specifications can be according to customers' requirements.
Customized Logo/brand printed.
Various colours available.

LED mouth mirror



1. Dimensions: 14 (L) x 1 (W) x 2.3 (H) cm
2. Rubber-in-plastic handle and switch button, easy to operation
3. Material: ABS
4. Use 3 piece of button battery
5. Mirror diameter: 23mm
Customized Logo/brand printed.
Various colours available.
10boxes/ctn
100pcs/bag 10bag/box 10boxes/ctn

Disposable calibration probe



Packing details: 50 ps/bag, 40 bag/ctn
G.W. 9,5 kg

USB Dental camera



Specifications

dental inspection camera
1. 1,480,000 pixels
2. 6 white intense LED
3. 1.8 m usb cable
4. with software

PRODUCTS / Poly sleeve


Barrier Film



Universal Barrier Film

Size:4"x6", 1200 sheets/roll with edge (5mm no adhesive-without edge (all with adhesive including edge)
Packing:1 roll/box, 12boxes/ctn
Dimensions: 40x37.5x36.5cm
G/N weight: 13kg/15.1kg/ctn
Color: clear, blue, pink

Dispenser



Product Type:
Oral Therapy Equipments & Accessories
Material: ABS
Certificate: CE,ISO13485,FDA
Application: Dental Areal
Future: Durable
Color: White

Plastic sleeve



Disposable dental sleeve ie ,dental chair cover ,full chair cover, half chair cover, ultrasonic scaler l, high&low speed handpiece , contra-angle handpiece , three way syringe ,x-ray tube, intraoral camera, tray cover,headrest,operating light handle etc

plastic X-ray Sensor Cover (Regular-thick 2.0s) 8" *1-5/8"	36boxes*500pcs/ctn	Plastic Tray Sleeves(Cover) (Extra-thick,2.5s) 10-1/2"*14"	6boxes*500pcs/ctn
Plastic X-ray Sensor Cover (Regular-thick 2.0s)10" *2-1/2"	36boxes*500pcs/ctn	Plastic X-Ray Sleeve (Regular-thick,2.0s) 30"*24"	250pcs/ctn
Plastic Scaler Sleeve(Regular-thick 2.0s) 1-3/5" * 8"	36boxes*500pcs/ctn	Plastic X-Ray Sleeve (Regular-thick,2.0s) 26"*15"	4boxes*250pcs/ctn
Disposable Intraoral Camera Sheath 9-4/5" * 2"	36boxes*100pcs/ctn	Plastic X-Ray Sleeve (Extra-thick,2.5s) 30"*24"	250pcs/ctn
Plastic Head Rest Cover (Regular-thick,2.0s) 11"*10"	12boxes*250pcs/ctn	Plastic X-Ray Sleeve (Extra-thick,2.5s) 26"*15"	4boxes*250pcs/ctn
Plastic Head Rest Cover (Regular-thick,2.0s) 14"*10"	12boxes*250pcs/ctn	Full Chair Sleeves (Regular-thick,2.0s) Full,29"*80"	3boxes*125pcs/ctn
Plastic Head Rest Cover (Extra-thick,2.5s) 11"*9-1/2"*2"	12boxes*250pcs/ctn	Full Chair Sleeves (Regular-thick,2.0s) Full,48"*56"	3boxes*150pcs/ctn
Plastic Head Rest Cover (Extra-thick,2.5s) 14"*9-1/2"*2"	12boxes*250pcs/ctn	Half Chair Sleeves (Regular-thick,2.0s) Half,27-1/2"*24"	3boxes*250pcs/ctn
Plastic Handpiece Sleeves (Regular-thick, 2.0s) 1"*8"(the tip length 4cm)for high-speed handpiece	36boxes*500pcs/ctn	Plastic Full Chair Sleeves(Extra-thick 2.5s) Full,29"*80"	3boxes*125pcs/ctn
Plastic Handpiece Sleeves (Regular-thick, 2.0s) 1-1/2"*8"(the tip length 4cm)for low speed, straight handpiece	36boxes*500pcs/ctn	Plastic Full Chair Sleeves(Extra-thick 2.5s) Full,48"*56"	3boxes*150pcs/ctn
Plastic Handpiece Sleeves (Regular-thick, 2.0s) 1-1/2"*10"(the tip length 4cm)	36boxes*500pcs/ctn	Plastic Half Chair Sleeves(Extra-thick 2.5s) Half,27-1/2"*24"	3boxes*250pcs/ctn
Plastic Handpiece Sleeves (Regular-thick, 2.0s) 1-1/2"*10"(the tip length 4cm) with slip at top	36boxes*500pcs/ctn	Plastic Full Chair Sleeves (thickness: 1.34s) Full,29"*80"	3boxes*125pcs/ctn
Plastic Handpiece Sleeves (Regular-thick, 2.0s) 1.3" * 8" (the tip length 4cm)	36boxes*500pcs/ctn	Plastic Full Chair Sleeves (thickness: 1.34s) Full,48"*56"	3boxes*150pcs/ctn
Plastic Handpiece Sleeves (Regular-thick, 2.0s) 7-3/10" * 5" (for low speed, curved handpiece)	36boxes*500pcs/ctn	Plastic Half Chair Sleeves (thickness: 1.34s) Half,27-1/2"*24"	3boxes*250pcs/ctn
Palstic Syringe Sleeves (Regular-thick, 2.0s) 2-1/2"*10"(transparent)	36boxes*500pcs/ctn	Plastic Light Sleeve for T-style Handle 4" * 5.75" (10.5cm X 15cm)w	36boxes*500pcs/ctn
Palstic Syringe Sleeves (Regular-thick, 2.0s) 2-1/2"*20")	36boxes*500pcs/ctn	Plastic Light Sleeve for T-style Handle 6"*2-9/10"	36boxes*500pcs/ctn
Palstic Syringe Sleeves (Regular-thick, 2.0s) 2-1/2"*10"(blue))	36boxes*500pcs/ctn	1 Full -Cover Barrier Film (10cm x 15cm x 1200 sheets) transparent,stick edges	8boxes*1rolls/ctn
Palstic Syringe Sleeves (Extra-thick,2.5s) 2-1/2"*10"	36boxes*500pcs/ctn	1 Full -Cover Barrier Film (10cm x 15cm x 1200 sheets) blue,stick edges	8boxes*1rolls/ctn
Palstic Syringe Sleeves (Extra-thick,2.5s) 2-1/2"*20"	36boxes*500pcs/ctn	1 Full -Cover Barrier Film (10cm x 15cm x 1200 sheets) transparent, non-stick edges	8boxes*1rolls/ctn
Palstic Syringe Sleeves (Extra-thick,2.5s) 2-1/2" *10" (hand tip cutting)	36boxes*500pcs/ctn	1Full -Cover Barrier Film (10cm x 15cm x 1200 shets) blue,non-stick edges	8boxes*1rolls/ctn
Palstic Syringe Sleeves (Extra-thick,2.5s) 2-1/2"*20"(hand tip cutting)	36boxes*500pcs/ctn	Full-cover Barrier Film Dispenser	8boxes*1 unit/ctn
Plastic Tray Sleeves(Cover) (Regular-thick,2.0s) 11-5/8"*16"	6boxes*500pcs/ctn	2 Air-motor &Syringes disposable shield (tubes)	10boxes*250pcs/ctn
Plastic Tray Sleeves(Cover) (Regular-thick,2.0s) 10-1/2"*14"	6boxes*500pcs/ctn	3 Tubing & turbines disposable shield (tubes)	10boxes*250pcs/ctn
Plastic Tray Sleeves(Cover) (Extra-thick,2.5s) 11-5/8"*16"	6boxes*500pcs/ctn	4 Headrests & X-ray heads disposable shield	8boxes*250pcs/ctn

PRODUCTS / Microbrush

Micro Applicator



1. Disposable, Non Absorbent, non-lint
Regular: Blue & light Blue (Dia:2.5mm);
Fine: Green & light Green (Dia:2.0mm); Ultra fine:
Purple & light Purple (Dia 1.5mm)
Cylinder: White & Black (Dia 1.2mm)
Package: 100pcs/bundle, 240 bundle / ctn

Microbrush divider.



Simply Micro Applicator
Size: 146*72*92
Packing details: 1 case/box, 50 boxes/ctn

Dowble microbrush



Double-ends brush
(pink, blue)
double-ends brush,
color: pink, blue 100pcs/
box

Microbrush tip with handle



Brush Applicator with handle
Handle, 100pcs/bag, 30bag/ctn
11-70032: Brush, 100pcs/bag

Microtip brush and handle



Microtip brush and handle
Color: Blue, Yellow, Green, Purple
2pcs handle (2pcs/box) + 400pcs
Microtip (100pcs/box)/kit,
70kits/ctn

Brush applicator




Bendable brush
Color: Blue, pink
Packing: 100pcs/drum, 4drum/
box, 32box/ctn

Brush applicator



1. Disposable
Description:
1. Chemical resistant fibers hold up to the
harsh dental materials
2. Bends close to tip for improved control and
ease of reach
3. Usage: whitening, fluoride, etchants, bonding
agents, sealants, desensitizers, adhesives /
resin. Color: Blue, red, green, yellow
Material: Medical Grade Plastic, Black Nylon
Brush Package: 144 Pcs/box

Adhesive Tip applicators



Usage
Pick up small item gently with
adhesive tip.
Place small item on prepared
surface. Twist to release small
item from adhesive.
Ideal for Veneers Inlay/Onlay
Cementation, Crowns, Bands,
Brackets, Clasps

PRODUCTS / Air-water tips

Air-water syringe



Air-water syringe
1pc air water syringe+2pcs tips
1 kit/box, 100box/ctn

Cover for air-water syringe tip



Air-water syringe tip cover

Color: white

Packing: 250pcs/bag, 40bag/ctn

Air-water syringe tip, metal core



3- water nozzle tips with metal core
Air-water syringe tip, metal core
white tube, metal core

Packing: 200pcs/box, 32box/ctn

Air-water syringe tip, full metal



Air-water syringe tip, full metal -3 way

1. Material: stainless steel.

2. Autoclavable

5pcs/bag

Air-water syringe tip, plastic colorfull tube



Air-water syringe tip, full plastic

colorful tube: blue, clear, white, red, purple, yellow, green are available

Packing: 250pcs/bag, 40bag/ctn

Air-water syringe tip, plastic colorfull core



Air-water syringe tip, full plastic with clear tube and colorful core:

blue, clear, white, red, purple, yellow, green are available

Packing: 250pcs/bag, 40bag/ctn

Air-water syringe tip, plastic, no core



white tube, no core
Packing: 250pcs/bag, 50bag/ctn

transparent tube, no core
Packing: 250pcs/bag, 50bag/ctn

Air-water syringe tip



straight, white
straight, assorted colours

pre-bent, white
pre-bent, assorted colours

PRODUCTS / Wedges, mixing tips


Fixing Plastic Wedges



Fixing Plastic Wedges with hole



Long handle wedges



Fixing Wooden Wedges




Mixing tips



Mixing Tip for Dental
High performance mixing tips
Different sizes/colors for different mixing needs.
Color-coding system precise & easy dispensing.

Amalgam carrier



Dental mixing intra oral tips



Dynamic Mixing tip



PRODUCTS / Dental laboratory

Flexible rubber mixing bowl



Mixing bowls: Flexible plastic material. Washable for repeated use
Usage: Flexible plastic material. Washable for repeated use.
Available of six different kinds of mixing bowls for different use.

	Size (Dia x h)	Volume	Color	Packing (ctn)	Weight (ctn)
L	130 * 96mm	690ml	Dark blue, black, clear pink, clear green	100pcs	13KGS
M	108 * 59mm	280ml	Dark blue, black, clear pink, clear green	180pcs	15KGS
S	84 * 50mm	150ml	dark blue, black, clear pink, clear green	250pcs	14KGS

Glass dappen dish



High quality glass dappen dish Autoclavable. Smooth surface for good work. Available Colour: Blue, White, Red, Green. This dappen dish is made of high quality, durable glass. Heavy and steady bottom to avoid spills. Suitable for working with small amounts of nail liquids, hair tints, esthetic glues and liquids. Perfect size for minimizing waste and contamination. Holds 5 ml

Size:	33*33mm
Color:	White, blue, green, amber
Packing:	12PCS/box, 18boxes/CTN
Carton size:	35*17*25.5cm
G.W.	10 KGS

No-stick mixing cups



Resistant to the high temperature of 134 C, Flexible silicone rubber materials. Five different size to fit your different jobs. Washable for repeated use. Surface cup bottom, with easy to put the bowls firmly on the smooth surface as glass, etc.

Disposable Dappen Dishes



Be put into using completely for time, may save time and be guarded against crossinfection. Be used for etching agent and binder.
Size: 23*32 mm
thermo stability
Colour: Pink, White, Blue, Yellow
Material: PP

Mixing Finger Dappen Dish



Convenient for mixing medicine, glue and so on
Perfect tools specially designed for Dental Lab work
Disposable, or compatible with the cold sterilization
Package Include:
200pcs in mixed colors

Rubbing bowl & rool



Material: Glass
Bowl: 48.3*28.1mm
Rool: 92.4*16mm
Package:
2set/box, 48boxes/ctn
410*270*225

Mixing spatula



Straight spatula, the handle is arabesquitic, easily of handling, flat and smooth fore-end of the spatula is better for mixing the medicine well. Autoclavable can be heated to 121
Specification: 195*32mm
Color: Pink, White, Blue, Yellow, Green

Mixing spatula



AutoClavable can be heated to 121, handle designing according to the shape of hands, closely and comfortable of holding. Big fore-end of the spatula, convenient for mixing up the medicine well. Small spatula, smooth surface, good effect of mixing well the medicine.
Specification: 185*33mm

PRODUCTS / Dental laboratory

Mixing Rod



Disposable mixing sticks, long handle, wide and flat tips, convenient of mixing well the medicine.
100pcs/bags, 70bags /ctn;
Carton size: 280*250*240

Disposable impression syringe



Curved Utility Syringe

Polypropylene plastic; self-enclosed medical piston. Ideal for teeth clearing, impression material injection, family care after operation, and the fluoride treatment.
Specification: 12ml
50pcs/bag, 25bag/ctn
500*360*430

Mixing Well



Disposable plastic mixing plates, 2 or 4 well, saving medicine, groove design on the two sides, convenient of handling.
180pcs/bags, 20bag/ ctn,
Carton size: 410*410*290mm
Specification: 50*50*10mm,
available 2 or 4 well

Impression Mixing Gun



Dental Impression Mixing Dispenser Dispensing Gun 50ml 1:1, 2:1 caulking Gun The dispensing gun is used to extrude impression materials and bite registration materials from cartridges. It consists of a release lever and plunger, which are coordinated simultaneously to load and remove the cartridge. Works with 1:1 and 10:1 cartridges.
Available in two styles for HP cartridges and old style cartridges. 50ml: 1:1, 10:1
Rational gear structure
Aircraft material applied, High hardness
Drop-back gear prevented
Package: Carton: L*W*H 50*40*35, 50pcs/ctn

Silicone impression material



Second step impression material for detail of supporting teeth (Injection type). Inlay, Onlay, Crown Bridge, Porcelain. Advantages. Excellent hydrophilic and perfect spread plus easy injection. Ideal flow guarantees perfect permeation into the sulcus without leakage. Best reproduction of details (25um) and lowest deformation. Vacuum mixing production system guarantees no air bubble.
Instructions Auto Type
1. Insert the cartridge into the dispenser.
2. Squeeze little bit until base and catalyst drain together without blocking.
3. Insert the mixing tip into the cartridge and turn it 90 to the clock wise direction.
4. Pull the trigger slowly for good mixing of the materials.

Endo Measuring Block plastic



Used for Root canal file measurement
Accurate settings. Durable construction. Fully autoclavable up to 135°C.
Size: L100mm*D24mm*H32mm
1pc/bag 1bag/box 300boxes/ctn
Color: Blue. Convenient plastic device for setting a precise working length. Plastic Economic alternative to autoclavable counterparts
Holds 32 measuring grades
Gauge from 12 to 27 mm.

Finger ring dental (endo ruler)



Product name: Measure scale
Function: for the measurement of root canal. 135C high temperature disinfection. Size: 35mm
Purpose: Pulp canal reamer measure. The span can be wear on the finger, it like a ring, the ring can be adjusted easily. Color: Blue (yellow purple can be selected)


Bur instrument box (bur holder)



bur Instruments Sterilizing Box
Fully autoclavable up to 135°C.
Color: yellow, blue, black
specification: 65*38*12mm
carton size:
410*270*235mm
weight 8.6kg

PRODUCTS / Dental products

Impression trays



High quality impression tray from Mediwish; strong, rigid distortion free plastic construction. Anatomically designed. Smooth edges for patient comfort. Excellent retention for impression materials whilst providing sufficient pressure for all details. 10 pcs/set. Big, Medium, Small, Side, front. Detail Order NO.: 1. large upper; 2. large lower; 3. medium upper; 4. medium lower; 5. small upper; 6. small lower; 7. upper right/lower left; 8. upper left/lower right; 9. anterior

Impressions trays steel



Disposable impression trays Strong, rid distortion. Anatomically designed. Smooth edges for patient comfort. Excellent retention for impression materials whilst providing sufficient pressure for all details

Autoclavable impression tray. Autoclavable up to 135 degree.

Dental Fiber Post



Diameter: 1.0mm, 1.2mm, 1.4mm, 1.6mm
Length: 17mm; Non-transparent, tapered; Radiopaque;
Fiber: 65%, Resin: 35%.

Absorbent Paper Points set



Gutta Percha Points and Absorbent Paper Points Mixed Round Box Packages Uniform and non-distorting crimp. Radiopacity. Long shelf life and freshness assured. Suitable for vertical and lateral condensation. Colour Coded Paper Points are made from specially formulated. Each of the Paper Point is individually hand jig rolled to ensure perfect size and firmness highly absorbent paper.

Size: 15#-40#

Dental X-ray film



Intra-Oral Dental Film ES-58.
150 count. Size: 31*41mm

X-ray film envelope



Size according with quotation

Gutta Percha Points




Description:
Hand rolled endodontic points used to fill a prepared root canal prior to tooth restoration.
Color Coded, ISO Size
All kinds of Accessory points.
Taper: 02
Size: 15-40, 45-80
Package: 120pcs/box , 10boxes/pack, 100box/carton

Tooth polisher



Material: synthetic rubber mixing with normal abrasive. Usage: For finishing and polishing titanium alloy/pure titanium/ nickel-chrome alloy/ cobalt- chrome alloy dental workpiece, such as for steel teeth/ inner crown/ steel frame/steel base finishing, smoothing and polishing. silicon polishers (special for porcelain/resin teeth). Material: silicon rubber mixing with normal abrasive. Usage: For smoothing and initial polishing porcelain/resin teeth.(super bright or without glaze polishing. suggest use RD2/0 and SD2/0 series product). silicon polisher (special for resin base). Material: silicon rubber mixing with normal abrasive. Usage: For smoothing and polishing the resin base. F series: for alloy workpieces. F series polishers are specially designed for alloy workpieces, used for various alloy tooth/alloy frames/alloy base plates grinding and polishing, also used for tooth crown & connecting points of the alloy workpieces fine finishing and polishing. G series: for gold & precious metal workpieces. G series polishers are specially designed for gold & precious metal workpieces, used for gold tooth grinding and polishing, also used for connecting points of the alloy workpieces fine finishing and polishing. H series: for ceramic workpieces. H series polishers are specially designed for ceramic workpieces, mainly for porcelain tooth surface/tooth crown/ tooth a lot/connecting points & edge of the workpieces fine finishing and polishing. K series: for resin workpieces. K series polishers are specially designed for resin workpieces, used for various resin tooth & hidden denture grinding and polishing; also used for tooth a lot/ connecting points of the workpieces/ various resin base plates/ crylic acid resin plates fine finishing and polishing. G series: for alloy & gold workpieces. G series polishers are specially designed for dental surgeries, mainly used for porcelain tooth surface/ tooth crown/ tooth a lot/ connecting points & edge of the workpieces fine finishing and polishing; also used for true tooth/ enamel/ resin workpieces finishing and polishing.



fine finishing and polishing. H series: for ceramic workpieces. H series polishers are specially designed for dental surgeries, mainly used for porcelain tooth surface/ tooth crown/ tooth a lot/ connecting points & edge of the workpieces fine finishing and polishing; also used for true tooth/ enamel/ resin workpieces finishing and polishing.

Silicon Rubber polishers Kits are specially designed for dental surgeries:

G series used for metal & precious metal workpieces. H series used for ceramic & resin workpieces as well as real tooth polishing.

PRODUCTS / Burs-Diamond & Carbide

Diamond Burs

Material:
Diamond grits mixing with plate metal binder

Usage: For finishing, carving and smoothing the slight parts of all kinds of dental workpiece. Such as steel toothless frame, inner steel crown (not suitable for the cut face of the crown), porcelain teeth/ resin teeth/ resin base.

Caution: Please use the products in allowed RPM. Or it may take the products damaged!

Order No.	Max RPM	Set B x L (mm)	ISO No.
514-12	25,000	1.2 x 1.2	514-12-01
524-12	25,000	1.2 x 1.2	524-12-01
534-12	25,000	1.2 x 1.2	534-12-01
544-12	25,000	1.2 x 1.2	544-12-01
514-16	25,000	1.6 x 1.6	514-16-01
524-16	25,000	1.6 x 1.6	524-16-01
534-16	25,000	1.6 x 1.6	534-16-01
544-16	25,000	1.6 x 1.6	544-16-01
514-20	25,000	2.0 x 2.0	514-20-01
524-20	25,000	2.0 x 2.0	524-20-01
534-20	25,000	2.0 x 2.0	534-20-01
544-20	25,000	2.0 x 2.0	544-20-01



Base on ISO Standards color and grit. Please refer quotation sheet bur's chart.

Red	514	Fine	F	50-60 um	#270-230
Blue	524	Standard	M	105-125 um	#140-120
Green	534	Coarse	C	125-150 um	#120-100
Black	544	Super Coarse	EC	150-180 um	#100-80

"CARBIDE BURS"

Packed 1 per unit

The same top quality carbide burs used in manufacturing CARBIDE BUR GOLDIES are also available uncoated as follows:



PRODUCTS / instrument containers

Bur holder (bur box) available in different styles (root canal dispenser)



Bur Holder 24 Slot



Bur Holder Box60 Slot Round



Bur Holder 24 Slot Round




Bur Magnetic Holder



Bur Holder 24 Slot



Bur Holder Box60 Slot Round



Bur Holder 24 Slot R



24 hole Bur Stand



PRODUCTS / instruments container

Endo boxes available in different styles



PRODUCTS / Propphy brush&cup

Propphy Brush junior bristle



Product Description:
disposable propphy brush
package: 100 pcs/ bag,
100 bags per carton

Propphy Brush flat nylon



Product Description:
disposable propphy brush
package: 100 pcs/ bag,
100 bags per carton

Propphy Brush tapered



Propphy Brush junior nylon



Propphy Angle




Angle propphy specifica-
tion:
Latex free, fit all standard
straight handpieces,
100% guarantee to last
an entire propphy, ribbed
and webbed propphy cups
Color: Soft cup(green),
Regular cup(blue), Hard
cup(purple)
Snap-on with serrated
knob, 100PCS/OPP bag
,10bags/ctn

Disposable Propphy Cups



Product Description:
with 4 angles or 6 anles
disposable propphy cup
Color: white, soft green, soft
blue, soft grey,
package: 100 pcs/ bag, 100
bags per carton



Tapered

Flat

Cup-shape

Dental Propphy Powder

Product Description:

High quality CE flavored
dental best propphy powder
best propphy powder

Three flavor For air polisher



Snap-on

Screw type

Latch type

PRODUCTS / Safety glasses

safety glasses



1.Material:Anti-scratch PC lens,Nylon frame
2.99.9% UV protection,anti-impact,anti-fog,anti-scratch
Color-black, blue, red,
Description- anti-fog/UV protection/anti-shock/anti-scratch/adjustable frame
Weight-33g/pcs, Package, one pc per polybag
Carton Size, 51*44*35CM, Weight-GW12KG NW11KG

safety glasses



1.Material:Anti-scratch PC lens,Nylon frame
2.99.9% UV protection,anti-impact,anti-fog,anti-scratch
Color-black, blue, red,
Description- anti-fog/UV protection/anti-shock/anti-scratch/adjustable frame
Weight-33g/pcs, Package, one pc per polybag
Carton Size, 51*44*35CM, Weight-GW12KG NW11KG

safety glasses



1.Material:Anti-scratch PC lens,Nylon frame
2.99.9% UV protection,anti-impact,anti-fog,anti-scratch
Color-black, blue, red,
Description- anti-fog/UV protection/anti-shock/anti-scratch/adjustable frame
Weight-33g/pcs, Package, one pc per polybag
Carton Size, 51*44*35CM, Weight-GW12KG NW11KG

safety glasses



1.Material:Anti-scratch PC lens,Nylon frame
2.99.9% UV protection,anti-impact,anti-fog,anti-scratch
Color-black, blue, red,
Description- anti-fog/UV protection/anti-shock/anti-scratch/adjustable frame
Weight-33g/pcs, Package, one pc per polybag
Carton Size, 51*44*35CM, Weight-GW12KG NW11KG

Shiled protective cover hood



Shiled protective cover hood
1.Material:Anti-scratch PC lens,Nylon frame 2.99.9% UV protection,anti-impact,anti-fog,anti-scratch
Color-transparent, anti-fog/UV protection/anti-shock/anti-scratch/adjustable frame, Weight 8.5kg/ctn ,Package-1shelf 10sheet, Carton Size-60*37*44cm, product size-sheet 277*210mm

Shiled protective cover hood



Shiled protective cover hood
1.Material:Anti-scratch PC lens,Nylon frame 2.99.9% UV protection,anti-impact,anti-fog,anti-scratch
Color-transparent, anti-fog/UV protection/anti-shock/anti-scratch/adjustable frame, Weight 8.5kg/ctn ,Package-1shelf 10sheet, Carton Size-60*37*44cm, product size-sheet 277*210mm

Shiled protective cover hood



Shiled protective cover hood
1.Material:Anti-scratch PC lens,Nylon frame 2.99.9% UV protection,anti-impact,anti-fog,anti-scratch
Color-transparent, anti-fog/UV protection/anti-shock/anti-scratch/adjustable frame, Weight 8.5kg/ctn ,Package-1shelf 10sheet, Cartor Size-60*37*44cm, product size-sheet 277*210mm

Shiled protective cover hood



Shiled protective cover hood
1.Material:Anti-scratch PC lens,Nylon



frame 2.99.9% UV protection,anti-impact,anti-fog,anti-scratch
Color-transparent, anti-fog/UV protection/anti-shock/anti-scratch/adjustable frame, Weight 8.5kg/ctn ,Package-1shelf 10sheet, Carton Size-60*37*44cm, product size-sheet 277*210mm

PRODUCTS / Check-retractor

rond



Expand mouth



fish tai



support



Unilatera




Lateral



Lateral



single span



double span



Lateral



support



support, autoclave



PRODUCTS / Ortodontic

Interdental floss pick
toothpick



Dental floss



Interdental brushes



Ortodontic Bracket



Dental stainless steel
orthodontic wire




Ligature tie



Toothbrushes



Elastic power chain



Dental air compressor



OIL FREE: without lubricated oil, harmless to the human body.
SUPER SILENT: Noise level lower than 56DB(A) to reduce noise pollution.
Power: 1100W(1.5HP)
Supply Voltage: 220V AC, 50HZ (220V AC, 60HZ 110V AC, 60HZ)
Air Displacement: 140L/min
Noise: 53db, Tank: 55L, Max Pressure: 8bar,
Weight: 65kgs
Dimension: 480*900*700mm

Dental Halogen Light Cure



PRODUCTS / Bandage, plasters, wound care products

Skin closure strip



Adhesive Wound Closures
Adhesive Wound Closures With An Exceptionally Porous Spun-Bond Fiber Backing, Includes Bilateral Nylon Filaments For High Tensile Strength. The Hypoallergenic Medical Grade Adhesive Affords Excellent Adhesion While Being Gentle To The Skin. Latex-Free. Wound Closure Strip(mm). Package MOQ(sheet) 3x75mm5pieces/sheet,1sheet/paper blister bag,1000sheets/carton5000; 6x75mm3pieces/sheet, 1sheet/paper blister bag, 1000sheets/carton 5000; 6x100mm 5pieces/sheet,1sheet/paper blister bag,1000sheets/carton 5000; 12x100mm 3pieces/sheet,1sheet/paper blister bag,1000sheets/carton 5000

Alginate dressing



Description:

CM Alginates is soft,highly absorbent dressing that quickly forms a hydrophilic gel to protect and cover the woundbed.

Outstanding absorbent capacity for fewer dressing changes,helps create moist wound healing environment,enhances patient comfort

1,Ca and Na exchange and turn to gel. the gel keep wound moist.

2,Stimulate blood platelet

agglutination,fast bleeding stop

Applications:

Pressure ulcers,arterial ulcers,venous ulcers,diabetic ulcers,superficial wounds such as cuts and abrasions,donor sites,post-operative wound,trauma wounds,other dermal lesions,help control minor bleeding.

Medical tape



Packing: 1pc/envelope ; 50pcs/box ; 20 boxes/ctn Material: non-woven or PU

Size: 5*7cm,8*10cm,10*10cm;15*10cm;20*10cm;25*10cm;30*10cm;35*10cm MOQ: 100000pcs

OEM: Provided, Size:Length(stretched): 4m,4.5m,5m, Width:4cm,5cm,6cm,7.5cm,8cm,10cm,15cm

Zinc Oxide Adhesive Plaster



Silk Plaster Features:

- 1) Silk material
- 2) Good adhesive ability
- 3) Easy to store, long storage life
- 4) Harmless to skin
- 5) Dimensions: 2.5x500cm/roll,5x500cm/roll,7.5x500cm/roll, other size available on requests, package: 12rolls/box, 30boxes/ctn

G.W./N.W.: 15/14KGS, carton size: 37x36.5x44.5cm
Container: 500ctns/20FT, 1100ctns/40HQ
Size:Length(stretched): 4m,4.5m,5m,
Width:4cm,5cm,6cm,7.5cm,8cm,10cm,15cm

PBT bandage, PBT Elastic Bandage



Bleached tabby PBT bandage
Material: 55% viscose: 45% cotton with woven edges
Weight g/sqm: 28g, 30g
Clip: With or without clips,elastic band clips or metal clips
Size: Length (stretched): 4m, 4.5m, 5m,
Width: 4cm, 5cm, 6cm, 7.5cm,8cm, 10cm, 15cm,
Basic packing:Individually packed in cellophane

Non-woven Cohesive Wrap Bandage



Cohesive elastic bandage with cotton and spandex woven material available
Adheres to itself, non-sticky to hair or skin;
No pins or clips needed

Porous, soft, allow skin to breath and comfortable to users; Protects primary dressings; Provide controlled compression;
Water resistant, will not loose by sweat or water; Easy hand tear, no scissors needed.

Color: available in, multiple colors to meet your application needs. Size: 2.5cmx5y,

5cmx5y, 7.5cmx5y, 10cmx5y, 15cmx15y available, also the others can be made according to your requirement. Fixing bandages for splints, monitors and IVs. Pressure bandages to promote circulation and healing. Compression bandages to help control swelling and stop bleeding. Industrial First Aid Bandages Size:Length(stretched): 4m,4.5m,5m, Width:4cm,5cm ,6cm,7.5cm,8cm,10cm,15cmBasic packing:Individually packed in cellophane.

Non-woven adhesive plasters



Packing: 1pc/envelope ; 50pcs/box ; 20 boxes/ctn

Material: non-woven or PU
Size: 5*7cm,8*10cm,10*10cm;15*10cm;20*10cm;25*10cm;30*10cm;35*10cm

MOQ: 100000pcs
OEM: Provided, Size:Length(stretched): 4m,4.5m,5m, Width:4cm,5cm ,6cm,7.5cm,8cm,10cm,15cm

Silk medical tape



Silk Plaster Features:

- 1) Silk material
- 2) Good adhesive ability
- 3) Easy to store, long storage life
- 4) Harmless to skin
- 5) Dimensions: 2.5x500cm/roll,5x500cm/roll,7.5x500cm/roll, other size available on requests, package: 12rolls/box, 30boxes/ctn

G.W./N.W.: 15/14KGS, carton size: 37x36.5x44.5cm
Container: 500ctns/20FT, 1100ctns/40HQ
Size:Length(stretched): 4m,4.5m,5m,
Width:4cm,5cm,6cm,7.5cm,8cm,10cm,15cm

PRODUCTS / Medical instrument

ID baby braselets



Identification tape



Nasal Cannula



Urine cup



Vaginal speculum



Sizes: Large, Medium, Small. (In compliance with international standrad)

Material: Polystyrene (short for: PS)
Packing: 1 PCS/bag, 100 bags/carton.
We have three types of bag for our customers choice. They are polyethylene bag (PE bag), polypropylene bag (PP bag), and paper-plastic bag. The carton size is: 52*39*32 cm
Function: They are widely used for vaginal examination and treatment.

First aid kit



71 pieces or 124 pieces or Customized contents

2x Triangula Bandage (96x96x136cm) ; 2x Disposable Glove (Large size); 2x Saline wound wash (15ml); 1x Plastic Tweezer (13cm)
2x PBT bandage (5cm x 450cm); 2x PBT bandage (7.5cm x 450cm)
2x PBT bandage (10cm x 450cm); 1x Non-woven Medical tape (2.5cm x 5m)
1x Scissor (15cm) ;1x Disposable CPR mask (One way valve)
6x Safty pin (3.7cm); 10x Sterile Nonwoven swab (7.5cm x 7.5cm)
10x Alcohol wipe (5cm x 5cm); 2x Trauma Pad (10cm x 20cm)
20x Band aid (1.9cm x 7.2cm); 4x Cleaning Wipe (6cm x 8cm)
1x Pencil; 1x Note Pad (A7); 1x Instruction manual (A4)
1x First Aid Kit Bag (24cm x 19cm x 8cm)

ENT kit



Item Name	Specifications	Q-ty
Tray	white, plastic	1 pc
Rinoscope(Naso-scope)	plastic	1 pc
Laryngeal mirror	plastic or ABS, size S	1 pc
Laryngeal mirror	plastic or ABS, size M	1pc
ENT curved forceps	galvanized	1 pc
Ear funnel speculum	black color, diameter 4.0 mm	1 pc
Ear funnel speculum	black color, diameter 4.6 mm	1 pc
Tongue Depressor	150°18°1.6mm	1 pc
Dental Bib	40°50cm	1 pcs
Ear Swab	wooden, 15cm or 10cm can be choosen	5pc

Packaging Details: Individual package into PE bag, 110sets per carton

Tourniquet Identification tape



Tweezers, forceps, sponge holder



1.Material: ABC , PP plastic, different sizes
meet customer ,s demand
13,5, 14,7, 19, 25 cm

Who we are

We are a manufacturer of medical goods, including medical packages, one-off dental goods and one-off surgical goods, with well-equipped testing equipment and strong technical force. With a wide range, good quality, reasonable prices and stylish designs, our products are extensively used in clinics, hospitals, beauty spas and other industries.

Quality produce



Fast delivery



All our products have ISO 13485 & CE Certificates

