

Mediwish® Tray Liners: Enhancing Sterilization

Mediwish® Tray Liners are designed to enhance sterilization processes by effectively absorbing excess moisture, preventing wet packs and ensuring the integrity of your instruments.

 by David Kang



Placement of Mediwish Tray Liners

Inside the Tray

Place liner directly under instruments. This absorbs moisture released during sterilization.

Underneath the Tray (Optional)

An additional liner under the tray adds extra protection. This prevents moisture from reaching the tray.

Both

Combining both methods offers maximum protection against moisture.



Key Benefits of Mediwish Tray Liners

- 1 Effective Moisture Absorption**
Quickly absorbs excess moisture.
- 2 Condensate Control**
Prevents condensate buildup, ensuring dry packs.
- 3 Instrument Protection**
Protects instruments from damage.
- 4 Sterility Assurance**
Creates a dry environment to prevent microbial growth.



Efficiency and Instrument Integrity

Efficiency

Faster drying times improve the efficiency of sterilization.

Instrument Integrity

Protection from moisture and damage ensures instrument longevity.

Variety of Mediwish Tray Liners



Various Sizes

Available to fit various tray configurations.



Materials

Cellulose, poly/cellulose, cellulose/crepe options.



Absorbency Levels

Different gsm options to suit needs.



Product Specifications

Product No.	Description	Dimensions	Quantity
MZS004A	Tray Liner, cellulose, 80gsm	200mm x 305mm (7.87" x 12")	3000
MZS005A	Tray Liner, cellulose, 80gsm	250mm x 250mm (9.84" x 9.84")	2000
MZS006A	Tray Liner, cellulose, 80gsm	250mm x 500mm (9.84" x 19.70")	1000
MZS007A	Tray Liner, cellulose, 80gsm	305mm x 410mm (12" x 16.14")	1000
MZS008A	Tray Liner, cellulose, 80gsm	457mm x 550mm (18" x 21.65")	500
MZS009A	Tray Liner, cellulose, 80gsm	305mm x 305mm (12" x 12")	2000
MZS010A	Tray Liner, cellulose, 80gsm	305mm x 457mm (12" x 18")	1000
MZS011A	Tray Liner, cellulose, 80gsm	305mm x 610mm (12" x 24")	1000
MZS012A	Tray Liner, cellulose, 80gsm	457mm x 457mm (18" x 18")	500
MZS013A	Tray Liner, cellulose, 80gsm	457mm x 610mm (18" x 24")	500
MZSP012018	Tray Liner, poly/cellulose, 45gsm	305mm x 457mm (12" x 18")	600
MZSP020025	Tray Liner, poly/cellulose, 45gsm	508mm x 635mm (20" x 25")	400
MZS202500	Tray Liner, cellulose/crepe, 107gsm (with folded sides)	350mm x 580mm (13.78" x 22.83")	500
MZS204500	Tray Liner, cellulose/crepe, 107gsm (with folded sides)	350mm x 350mm (13.78" x 13.78")	500
MZS203500	Tray Liner, cellulose/crepe, 107gsm (with folded sides)	400mm x 600mm (15.75" x 23.62")	500





Material Composition and Safety

1 Durable Cellulose

High-quality cellulose material for effective absorption.

2 Low Chloride and Sulfate

Manufactured with low chloride and sulfate content for safety.

3 Free from Optical Dyes

Ensures compatibility with delicate instruments and processes.

Conclusion: Choosing Mediwish Tray Liners

Mediwish Tray Liners enhance sterilization. They improve efficiency, protect instruments, and ensure sterility. Choose the right liner for your needs. Contact your representative to learn more.

