

Facility	Department	Cleaner ID	Make / Model	Tank volume	Frequency	Detergent / conc.	Reference position

Defined routine method: empty basket/tray; bath prepared and degassed; 5-minute uninterrupted cycle; read immediately (preferably within 1 minute, no later than 5 minutes).

No.	Date	Cleaner ID	REF / LOT / EXP	Vial / Pos.	Temp °C	Cycle	Time to yellow	Result	Operator	Corrective action / sign-off
1								A / U / F		
2								A / U / F		
3								A / U / F		
4								A / U / F		
5								A / U / F		
6								A / U / F		
7								A / U / F		
8								A / U / F		
9								A / U / F		
10								A / U / F		
11								A / U / F		
12								A / U / F		
13								A / U / F		
14								A / U / F		
15								A / U / F		
16								A / U / F		
17								A / U / F		
18								A / U / F		
19								A / U / F		
20								A / U / F		
21								A / U / F		
22								A / U / F		
23								A / U / F		
24								A / U / F		
25								A / U / F		

Result codes: A = Acceptable yellow by 5 min; U = Unsatisfactory green/light green at 5 min; F = Fail / no visible change.

Yellow only after >5 to 7 min: record DELAYED in corrective-action field; investigate and repeat.

Any later color drift does not replace the result recorded at the prescribed read time.